



Recombinant Phosphoribosyl-AMP cyclohydrolase (hisI)

Product Code	CSB-MP700543MEZ
Abbreviation	hisI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q50837
Product Type	Recombinant Protein
Immunogen Species	Methanococcus vanniellii
Purity	>85% (SDS-PAGE)
Sequence	MGIKDIDIKE NFGKIVQNMD LKFRKIDDKE LLIAIAIDKY KNVLMTAFMD KESLKMTLKT GLMHYFSTSR NKIWMKGEES KNVQKVLEVF KDCDGDALLF IVEQTGWACH EGYMSCFHNK VDLNTGNSTV IGDKLD
Source	Mammalian cell
Target Names	hisI
Protein Names	Recommended name: Phosphoribosyl-AMP cyclohydrolase Short name= PRA-CH EC= 3.5.4.19
Expression Region	1-136
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.