



Recombinant *Oryza sativa* subsp. *japonica* Tubby-like F-box protein 13 (TULP13)

Product Code	CSB-YP700658OFG
Abbreviation	TULP13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53PP5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSFRSIVRDV RDGFGSLR R GFEVRILGHR RGKSHGAVHE LHDPVPVIQS SCWASLPPEL LRDIERLEE SEATWPSRKH VVACAGVCRT WREMCKEIVK NPELCGKITF PISLRQPGPR DGTMQCFIRR DKSTQTYLY LSLGSAVLVD NGKFLLSAKR NWHATCTEYV ISMNANNLSR STNTNIGKLR SNFLGTFKVI YDHTPYNAT SDSQSGKTSR RFSNKGTA KH PCSTYSIANI SYELNVFGTR GPRRMCCLMH SIPASSLEAG GTVPSQPDSI LAHSLNESSF RSVSFSKSSV MDHSMHFSSA QFSDISIGDG PRIGGRVLS D DEECKETPLI LQNKAPRWHE QLQCWCLNFR GRVTVASVKN FQLIAATQPA AGAPTPSQPV PPPPEHDKV ILQFGKVAKD MFTMDYHYPL SAFQAF AISL SSFDTKLACE
Source	Yeast
Target Names	TULP13
Protein Names	Recommended name: Tubby-like F-box protein 13 Short name= OsTLP13 Alternative name(s): Tubby-like F-box protein 6 Short name= OsTLP6
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.