



Recombinant Dictyostelium discoideum Casein kinase II subunit beta (csnk2b)

Product Code	CSB-EP700834DKK-B
Abbreviation	csnk2b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54PF1
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MDRKYKDTSS EEEDSLNEDE MAWIPWYCNL KGNEFFATID EDYINDDFNL TGLSSLVQFY DSALGIILDS DPDDPLSEDQ QEALERSADI LYGLIHARYI LTPKGLAHMH EKFKKAEFGR CPRVFCQNP VLPVGLADMV GVDTVKVCYCP RCNDIFNPKY RRHSHIDGAY FGTTFPHELL ITYPELIPTK PPQQYIPKIY GFKIHSARE RQLQIQKKN SLSNNNQNNQ NNNNNNNNNN NNNNNNNNNN NNNNQQNNNN QQNNNTNK
Source	E.coli
Target Names	csnk2b
Protein Names	Recommended name: Casein kinase II subunit beta Short name= CK II beta
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.