



Recombinant Dictyostelium discoideum Ubiquitin fusion degradation protein 1 homolog (ufd1)

Product Code	CSB-MP701013DKK
Abbreviation	ufd1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55BK0
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MNLDEYVRSL RGEVMGHHHG GPGRYEQKFK AFPISFLPKE KHSLESGGKI LLPPSALNAL SRLNIQYPML FEISNPISGK KSHCGVLEFI AEEGICYLPL WMMQNLQLKE GEFIDIKNAT LAKGTFVKIQ PRTSNFIDIS NPKAVLENSL RKFATLTKDD EIMIDYNNTK YYLKVVLPK ANAISIEAD VSVDFAPPMD SKEATSPSTS SPGSHVSGPS KGLTFGPAST SAKPIPGGKK KKDESDSDSD SDDEPKFKAF AGTGARLDGK VGTPLGTSPK TLNTNNNNNN NNNNNNNNNN NNNNNSNNNN TKNEDTKFKS FSGKGHSLKD
Source	Mammalian cell
Target Names	ufd1
Protein Names	Recommended name: Ubiquitin fusion degradation protein 1 homolog
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.