



Recombinant Danio rerio Probable arylformamidase (afmid)

Product Code	CSB-YP701098DIL
Abbreviation	afmid
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q566U4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSDWKSMSKD ELELQYSPSR WSPRMSADDV IRAHVSALKS GTERARSVTQ TLMDVPYGEA EGEKLDVYLP SSSSPDVPLV IYFHGGYWQF LSKDESGFLA VPLVQKGAVV VAVGYSIAPK GDMDLMVSQV RRSVVSVIQQ YSHISGLYLC GHSAGAHLAA MVLSTDWTQY DVSPKIKGAF LVSGIYDLQP ILSTYVNEPL KMTEEVALRN SPSRFLSQLK VSSASCDITV AVAQNDSPEF RKQSEEYYKA LESAGLKVSF EDVSETDHFS IIEQLVDENY LLTKLLLKMI GKS
Source	Yeast
Target Names	afmid
Protein Names	Recommended name: Probable arylformamidase EC= 3.5.1.9 Alternative name(s): Kynurenine formamidase Short name= KF
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.