



Recombinant Human Coiled-coil domain-containing protein 157 (CCDC157)

Product Code	CSB-EP701113HU-B
Abbreviation	CCDC157
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q569K6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAHLLGSQAC MESLRTDLTD LQGAIVDVFS RAGPVRFPSW KFPDRMACDL DMVALLEHYD HVPGDPEFTQ LSHAVLLELV IDRLLLLLQS CMSYLENLGS EQMMPPAQAA GPCMSVGLTV RRFWDSLLRL GTLHQQPLPQ KGANQRETPT SKPTTKGEP A RSPEYLTTKL IKPSSPVLGL PQTCEPESI PVRASLQFPA TTFKNTRSVH SQTLETALVP CDACASVQGS LQKVGKVVIS LCQSQNLPSS LGQFQQLVQD SMGLRPLPAA TVGRWAAEQR KDLTRLSKHV EALRAQLEEA EGQKDGLRKQ AGKLEQALKQ EQGARRRQAE EDEQCLSEWE HDKQQLLTET SDLKTKMATL ERELKQQRRES TQAVEAKAQQ LQEEGERRAA AERQVQQL EE QVQQLEAQVQ LLVGRLEGAG QQVCWASTEL DKEKARVDSM VRHQESLQAK QRALLKQLDS LDQEREELRG SLDEAEAQRA RVEEQLQSER EQGQCQLRAQ QELLQSLQRE KQGLEQATTD LRLTILELER ELEELKERER LLVAFPD LHR PTETQIHGGR SSSVESQITC PTDSGNVTDH MERQVQSNDI RIRVLQEENG RLQSMLSKIR EVAQQGGLKL IPQDRLWSPS SKGTQGATPP VQAKSTSPGP LGRQHLPSSR TGRTLLGQPC TSPPRQPCTS PPRQPCTSP RQPCTSPSRQ PCSQPSKSL L EGVTHLDTCT QNPIKVLVRL RKRLSPGRGQ ASSAHQPQER PM
Source	E.coli
Target Names	CCDC157
Protein Names	Recommended name: Coiled-coil domain-containing protein 157
Expression Region	1-752
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.