



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 6 (At1g19860)

Product Code	CSB-MP701130DOA
Abbreviation	At1g19860
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q56XU4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MRALHKRVSW PPDFKLCQVR LFISESPSQ VGSESQDHLQ AKSPLASHPS DDNLPPGGFGG PFSANESQIK LSDIPLIKWK CSVQILLDRE WKVVAGDESK EVEAQNEREL RVLEAFYPGA SSIPPNPSVP ADVEDSHHDD QQTIVIPILP VEDDDIAMDS ASDFPTQSGV DVGTEPSITD ENTSTSSTLP AGPDIMAALS AISNSKEQGS MIDQDLLIKI LSNPKLVENL VANRGSAGSV SSNTSSLYSS STHEANGVVT TAPISSNGQF YAQPPITHIP PMAYTPHAPQ DQPNYGAPPA RDASYYNLI QQHGGDRQET PPVQHLGYRY NLQPGGGPNP EMVNSSNNNQ RPRDSKPKIM KACMYFNSAR GCRHGANCY QHDATPYQPR NLNNGNINTS DMQNAKRMRF DRD
Source	Mammalian cell
Target Names	At1g19860
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 6 Short name= AtC3H6
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.