



# Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 6 (At1g19860)

<b>Product Code</b>	CSB-YP701130DOA
<b>Abbreviation</b>	At1g19860
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q56XU4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRALHKRVSW PPDFKLCQVR LFISEDSQSQ VGSESQDHLQ AKSPLASHPS DDNLPPGGFGG PFSANESQIK LSDIPLIKWK CSVQILLDRE WKVVAGDESK EVEAQNEREL RVLEAFYPGA SSIPPNPSVP ADVEDSHHDD QQTIVIPILP VEDDDIAMDS ASDFPTQSGV DVGTEPSITD ENTSTSSTLP AGPDIMAALS AISNSKEQGS MIDQDLLIKI LSNPKLVENL VANRGSAGSV SSNTSSLYSS STHEANGVVT TAPISSNGQF YAQPPITHIP PMAYTPHAPQ DQPNYGAPPA RDASYYNLI QQHGGDRQET PPVQHLGYRY NLQPGGGPNP EMVNSSNNQ RPRDSKPKIM KACMYFNSAR GCRHGANCY QHDATPYQPR NLNNGNINTS DMQNAKRMRF DRD
<b>Source</b>	Yeast
<b>Target Names</b>	At1g19860
<b>Protein Names</b>	Recommended name: Zinc finger CCCH domain-containing protein 6 Short name= AtC3H6
<b>Expression Region</b>	1-413
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.