



# Recombinant Human Myeloid zinc finger 1 (MZF1)

<b>Product Code</b>	CSB-BP015389HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P28698
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MRPAVLGSPD RAPPEDEGPV MVKLEDSEEE GEAALWDPGP EAARLRFRCF  RYEEATGPQE ALAQLRELCR QWLRPEVRSK EQMLELLVLE QFLGALPPEI  QARVQGQRPQ SPEEAAALVD GLRREPGGPR RWVTVQVQQQ  EVLSEKMEPS SFQPLPETEP PTPEPGPKTP PRTMQESPLG LQVKEESEVT  EDSDFLESGP LAATQESVPT LLPEEAQRCG TVLDQIFPHS KTGPEGPSWR  EHPRALWHEE AGGIFSPGFA LQLGSISAGP GSVSPHLHVP WDLGMAGLSG  QIQSPSREGG FAHALLLPSD LRSEQDPTDE DPCRGVGPAL ITTRWRSRPG  RSRGRPSTGG GVVRRGRCDV CGKVFSQRSN LLRHQKIHTG  ERPFVCSECG RSFSRSSHLL RHQLTHTEER PFVCGDCGQG FVRSARLEEH  RRVHTGEQPF RCAECGQSFR QRSNLLQHQR IHGDPPGPGA  KPPAPPGAPE PPGPFPCSEC RESFARRAVL LEHQAVHTGD KSFGCVECCE  RFGRRSVLLQ HRRVHSGERP FACAECGQSF RQRSNLTQHR RIHTGERPFA  CAECGKAFRQ RPTLTQHRLV HTGEKPFACP ECGQRFQSRL KLTRHQRTHT  GEKPYHCGEC GLGFTQVSRL TEHQRIHTGE RPFACPECGQ SFRQHANLTQ  HRRIHTGERP YACPECGKAF RQRPTLTQHL RTHRREKPFACQDCGRRFHQ  STKLIQHQRV HSAE</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	MZF1
<b>Protein Names</b>	Recommended name: Myeloid zinc finger 1 Short name= MZF-1 Alternative name(s): Zinc finger and SCAN domain-containing protein 6 Zinc finger protein 42
<b>Expression Region</b>	1-734
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.