



Recombinant Human Myeloid zinc finger 1 (MZF1)

Product Code	CSB-MP015389HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28698
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MRPAVLGSPD RAPPEDEGPV MVKLEDSEEE GEAALWDPGP EAARLRFRCF RYEEATGPQE ALAQLRELCR QWLRPEVRSK EQMLELLVLE QFLGALPPEI QARVQGQRPQ SPEEAAALVD GLRREPGGPR RWVTVQVQQQ EVLSEKMEPS SFQPLPETEP PTPEPGPKTP PRMQESPLG LQVKEESEVT EDSDFLESGP LAATQESVPT LLPEEAQRCG TVLDQIFPHS KTGPEGPSWR EHPRALWHEE AGGIFSPGFA LQLGSISAGP GSVSPHLHVP WDLGMAGLSG QIQSPSREGG FAHALLLPSD LRSEQDPTDE DPCRGVGPAL ITTRWRSRPG RSRGRPSTGG GVVRRGRCDV CGKVFSQRSN LLRHQKIHTG ERPFVCSECG RSFSRSSHLL RHQLTHTEER PFVCGDCGQG FVRSARLEEH RRVHTGEQPF RCAECGQSFR QRSNLLQHQR IHGDPPGPGA KPPAPPGAPE PPGPFPCSEC RESFARRAVL LEHQAVHTGD KSFGCVECCE RFGRRSVLLQ HRRVHSGERP FACAECGQSF RQRSNLQHR RIHTGERPFA CAECGKAFRQ RPTLTQHHRV HTGEKPFACP ECGQRFSQRL KLTRHQRTHT GEKPYHCGEC GLGFTQVSRL TEHQRIHTGE RPFACPECGQ SFRQHANLTQ HRRIHTGERP YACPECGKAF RQRPTLTQHL RTHRREKPFACQDCGRRFHQ STKLIQHQRV HSAE</p>
Source	Mammalian cell
Target Names	MZF1
Protein Names	Recommended name: Myeloid zinc finger 1 Short name= MZF-1 Alternative name(s): Zinc finger and SCAN domain-containing protein 6 Zinc finger protein 42
Expression Region	1-734
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.