



Recombinant Human F-box only protein 48 (FBXO48)

Product Code	CSB-BP702041HU
Abbreviation	FBXO48
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5FWF7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MHKNSKRNNN LRVSHTEANS VDAEKEKNES QNNFFELLPA EITFKIFSQL DIRSLCRASL TCRSWNDTIR NSDSLWKPHC MTVRAVCRRE IDDDLESGYS WRVILLRNYQ KSKVKHEWLS GRYSNICSPI SLPEKIMYPM DADTWGEILE AELER
Source	Baculovirus
Target Names	FBXO48
Protein Names	Recommended name: F-box only protein 48
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.