



# Recombinant Human Transcription factor Dp family member 3 (TFDP3)

<b>Product Code</b>	CSB-EP702134HU-B
<b>Abbreviation</b>	TFDP3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5H9I0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAKYVSLTEA NEELKVLMD E NQTSRPVAVH TSTVNPLGKQ LLPKTFGQSS VNIDQQVVIG MPQRPAASNI PVVGSPNPPS THFASQNHQS YSSPPWAGQH NRKGEKNGMG LCRLSMKVWE TVQRKGTTC QEVV GELVAK FRAASNHASP NESAYDVKNI KRRTYDALNV LMAMNIISRE KKKIKWIGLT TNSAQNCQNL RVERQKRLER IKQKQSELQQ LILQQIAFKN LVLRNQYVEE QVSQRPLPNS VIHVPFIIS SSKKTVINCS ISDDKSEYLF KFSSFEIHD DTEVLMWMMGM TFGLESGSCS AEDLKMARNL VPKALEPYVT EMAQGTGGV FTTAGRSNG TWLSASDLTN IAIGMLATSS GGSQYSGSRV ETPAVEEEEE EDNNDLSE NDEDD
<b>Source</b>	E.coli
<b>Target Names</b>	TFDP3
<b>Protein Names</b>	Recommended name: Transcription factor Dp family member 3 Alternative name(s): Cancer/testis antigen 30 Short name= CT30 Hepatocellular carcinoma-associated antigen 661
<b>Expression Region</b>	1-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.