



Recombinant Human Protein FAM76B (FAM76B)

Product Code	CSB-EP702405HU-B
Abbreviation	FAM76B
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5HYJ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AASALYACT KCTQRYPFEE LSQGQQLCKE CRIAHPIVKC TYCRSEFQQE SKTNTICKKC AQNVKQFGTP KPCQYCNIIA AFIGTKCQRC TNSEKKYGGP QTCEQCKQQC AFDRKEEGRR KVDGKLLCWL CTLSYKRVLQ KTKEQRKSLG SSHSNSSSSSS LTEKDQHHPK HHHHHHHHHH RHSSSHHKIS NLSPEEEQGL WKQSHKSSAT IQNETPKKKP KLESKPSNGD SSSINQSADS GGTDNFVLIS QLKEEVMSLK RLLQQRDQTI LEKDKKLTEL KADFQYQESN LR TKMNSMEK AHKETVEQLQ AKNRELLKQV AALSKGKKFD KSGSILTSP
Source	E.coli
Target Names	FAM76B
Protein Names	Recommended name: Protein FAM76B
Expression Region	2-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>