



Recombinant Human NHL repeat-containing protein 3 (NHLRC3)

Product Code	CSB-YP702496HU
Abbreviation	NHLRC3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JS37
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GSPVL RNFTFAVSWR TEKILYRLDV GWPKEHPEYFT GTTFCVAVDS LNGLVYIGQR GDNIPKILVF TEDGYFLRAW NYTVDPHGI FFASTLYEQS VWITDVGSGF FGHTVKKYSS FGDLVQVLGT PGKKGTS LNP LQFDNPAELY VEDTGDIYIV DGDGGLNNRL IKLSQDFMIL WLHGNGTGP AKFNIPHSVT LDSAGRVWVA DRGNKRIQVF DKDTGEWLGA WNNCFTEEGP SSVRFTPDGK YLIVAQLNLS RLSVVAAPPV GSIGECSVIS TIQLADQVLP HLLVDKRTG AVYVAEIGAK QVQKYVPLNS YVPSFGS
Source	Yeast
Target Names	NHLRC3
Protein Names	Recommended name: NHL repeat-containing protein 3
Expression Region	26-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.