



Recombinant Pig Perilipin-2 (PLIN2)

Product Code	CSB-MP703003PI
Abbreviation	PLIN2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4PLW0
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	ASVAVESQP SVVTRVANLP LVSSTYDLVS SAYISTKDQY PYLKSLCEMA EKGVKTITSV AMSGALPIIQ KLEPQIAIAN TYACKGLDRI EEKLPILNQP TNQVVANAKG AVTGAKDAMT TTVTGAkdCV ASTITEVVDK TKEAVTGSVE KTKSVVNGSI NTVLGGRMMQ LVSSGVEKAF TKSELLVDQY LPLTEEELEK EAKKVEGFDM VQKPSYYIRL GSLSTKLRSR AYQQALTRVK EVKQKSQETI SQLHSTVNLI EFARKNVHNA NQKIQGTQDK LYLSWVEWKR SIGYDDTDES HCAEHIESRT LAIARNLTQQ LQTTCHTLVS NIQGLPQNIH DQANHLGVMA GDIYSVFHNA SSFKEMSDGL LSSSKGQLQK MKESLDDVMD YLVNNTPLNW LVGPFYPQLT ESQDAQSRGA ENTTPETQQ PETKRIKPAP ASSAWGSQSG DTSCTVATC
Source	Mammalian cell
Target Names	PLIN2
Protein Names	Recommended name: Perilipin-2 Alternative name(s): Adipophilin Adipose differentiation-related protein Short name= ADRP
Expression Region	2-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.