



Recombinant *Xenopus laevis* Polycomb group RING finger protein 1 (pcgf1)

Product Code	CSB-BP703073XBE
Abbreviation	pcgf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4QR06
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MASQGPLVIA MRLRNQLQSV YKMDPLRNEE VVKV KIKELN EHIVCYLCAG YFIDATTITE CLHTFCKSCI VKYLQTSKYC PLCNIKI HET QPLLNLK LDR VMQDIVYKLV PGLQENEDSR IRDFYHSRGL ERVLQPSAVE DSVG DVSQLS LSLAVSQKTS HYRND EHV C LCLEKVSSGK DKKK FILQQK YVRC SVRSEI RHLRRVLSHR LSAPLAQVQL LIDNKVLPDH MTMKQLWLMH WYGK PAPLVL LYSVKEKRR
Source	Baculovirus
Target Names	pcgf1
Protein Names	Recommended name: Polycomb group RING finger protein 1
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.