



Recombinant Danio rerio Breast cancer metastasis-suppressor 1-like protein-A (brms1la)

Product Code	CSB-MP703295DIL
Abbreviation	brms1la
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4V8V1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MPVHSREKKE SNHEEMEVDV AEQEGSSSED EDTESSSVSE DGESSEMDDE DCERRRMECI DEMTNLEKQF TDLKDQLYKE RLNQVDAKLQ EVMSGKAAEY LEPLATLQEN MQIRTKVAGI YRELCLESVK NKYDCETQAA FQHWSEKLL LFDTVQTELE EKIRRLEEDR HSIDITSELW NDELQSRKNK KKDPFSPDKK KKPVVVSGPY IVYMLQQLDI LEDWTAIRKA MATLGPHRVK TDASSKIDKL QHNARSEDGR LFYDGEWYSR GQAITIDKKD EYPTSAVITT INHDEVWFKR VDGSKSKLYV SQLQKGKQTV KHA
Source	Mammalian cell
Target Names	brms1la
Protein Names	Recommended name: Breast cancer metastasis-suppressor 1-like protein-A
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.