



# Recombinant Mouse Interleukin-28A (IL28a)

<b>Product Code</b>	CSB-EP703302MO-B
<b>Abbreviation</b>	IL28a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q4VK74
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	D PVPRATRLPV EAKDCHIAQF KSLSPKELQA FKKAKDAIEK RLLEKDLRCS SHLFPRAWDL KQLQVQERPK ALQAEVALTL KVWENMTDSA LATILGQPLH TLSHIHSQLQ TCTQLQATAE PRSPSRRLSR WLHRLQEAQS KETPGCLEAS VTSNLFRLLT RDLKCVANGD QCV
<b>Source</b>	E.coli
<b>Target Names</b>	Ifnl2
<b>Protein Names</b>	Recommended name: Interleukin-28A Short name= IL-28A Alternative name(s): Interferon lambda-2 Short name= IFN-lambda-2
<b>Expression Region</b>	20-193
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a cytokine distantly related to type I interferons and the IL-10 family. This gene, interleukin 28B (IL28B), and interleukin 29 (IL29) are three closely related cytokine genes that form a cytokine gene cluster on a chromosomal region mapped to 19q13. Expression of the cytokines encoded by the three genes can be induced by viral infection. All three cytokines have been shown to interact with a heterodimeric class II cytokine receptor that consists of interleukin 10 receptor, beta (IL10RB) and interleukin 28 receptor, alpha (IL28RA).
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.