



# Recombinant *Xenopus tropicalis* Phosphatidylinositol 4-kinase type 2-alpha (pi4k2a)

<b>Product Code</b>	CSB-BP703498XBF
<b>Abbreviation</b>	pi4k2a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q505I0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDETSPLVSP DRDQTDYSYQ SQCSPGVPVC LSPNGRFSAV PGVVVRIPGS ATSPIRGSAA SGPSPPGSPC DRERQPLLER SQTRGAAAQA ERERNKFPDD PEFAEUVKKA EKAIVRDILP ERISQGSSGS YFVKNEQGEI IAVFKPKNEE PYGQLNPKWT KWLQKLCCPC CFGRDCLVLN QGYLSEAGAS LVDQKLELNI VPRTKVVFLA SETFNYS AID RVKSRGKRLA LEKVPKVGQR FNRIGLPPKV GSFQLFVKGY KDADYWLRRF EADPLPENTN RQLQLQFERL VVLDYIIRNT DRGNDNWLK YDCPMDSASA RDDWVMVKEP VIKIAAIDNG LAFPLKHPDS WRAYPFYWAW LPQAKIQFSQ EIKDLILPKI SDPNFVKDLE EDLYELFKRD PGFDRGQFRK QIAVMRGQIL NLQALKDGGK SPLQLVQTPP VIVETARSHQ KSTSESYTQS FQSRKPFFSW W
<b>Source</b>	Baculovirus
<b>Target Names</b>	pi4k2a
<b>Protein Names</b>	Recommended name: Phosphatidylinositol 4-kinase type 2-alpha EC= 2.7.1.67 Alternative name(s): Phosphatidylinositol 4-kinase type II-alpha
<b>Expression Region</b>	1-471
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.