



Recombinant Human Tigger transposable element-derived protein 5 (TIGD5)

Product Code	CSB-MP703611HU
Abbreviation	TIGD5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53EQ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MYPAGPPAGP VPRRGRRPLP GPPAPAPAPV PAARPPPPAP GPRPRVAVKM AFRKAYSIKD KLQAIERVKG GERQASVCRD FGVPGGTLRG WLKDEPKLRW FLEQLGGEVG TQRKKMRLAN EEEIDRAVYA WFLALRQHGV PLSGPLIQAQ AEFARQIYG PECTFKASHG WFWRWQKRHG ISSQRFYGEA GPPAPSPAPG PPVKEEPALP SGAGPLPDRA PAPPPPAEGG YGDEQIYSAS VTGLYWKLLP EQAAPPGAGD PGAGGCGRRW RGDRVTVLLA ANLTGSHKLG PLVIGRLPDP PSLRHHNQDK FPASYRYSPD AWLSRPLLRG WFFEEFVPGV KRYLRRSCLQ QKAVLLVAHP PCPSPAASMP ALDSEDAPVR CRPEPLGPPE ELQTPDGAVR VLFLSKGSSR AHIPAPLEQG VVAAFKQLYK RELLRLAVSC ASGSPLDFMR SFMLKDMLYL AGLSWDLVQA GSIERCWLLG LRAAFEPAPG EDSAGQPAQA EEAAEHSRVL SDLTHLAALA YKCLAPEEVA EWLHLDDDDGG PPEGCREEVG PALPPAAPPA PASLPSAMGG GEDEEEATDY GGTSVPTAGE AVRGLEALR WLENQDPREV GPLRLVQLRS LISMARRLGG IGHTPAGPYD GV
Source	Mammalian cell
Target Names	TIGD5
Protein Names	Recommended name: Tigger transposable element-derived protein 5
Expression Region	1-642
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.