



Recombinant *Oryza sativa* subsp. *japonica* Thaumatococcus-like protein (Os11g0706600, LOC_Os11g47944)

Product Code	CSB-EP703617OFG
Abbreviation	Os11g0706600, LOC_Os11g47944
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53MB8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	C GGGGGNGGIQ LIMVNCGES VWPGLLGTAG HPTPQSGGFH LGAGEEAALE VPAGWSGRVW PRRGCSFDSR GRGSCATGDC GGVLRNCGAA GATPATVVEM TLGTSASAMH FYDVSLVDGF NAPVSMAAVG GGVGCATAAC GADVNVCCPS ALEVRDREGR VAGCRSACRA MGGDYRCCTG DYASPSACRP TIFSHLFKAI CPRAYSAYD DATSLNRCHA KRYLITFCPP QPS
Source	<i>E. coli</i>
Target Names	Os11g0706600
Protein Names	Recommended name: Thaumatococcus-like protein
Expression Region	20-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.