



Recombinant Human HCLS1-binding protein 3 (HS1BP3)

Product Code	CSB-MP703623HU
Abbreviation	HS1BP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53T59
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQSPAVLVTS RRLQNAHTGL DLTVPQHQEV RGKMMSGHVE YQILVVTRLA AFKSAKHRPE DVVQFLVSKK YSEIEEFYQK LSSRYAAASL PPLPRKVLV FV GESDIRERRA VFNEILRCVS KDAELAGSPE LLEFLGTRSP GAAGLTSRDS SVLDGTDSQT GNDEEAFDFF EEQDQVAEEG PPVQSLKGED AEESLEEEEA LDPLGIMRSK KPKKHPKVAV KAKPSPRLTI FDEEVPDEG LFGPGRKLSP QDPSEDVSSV DPLKLFDDPD LGGAIPLGDS LLLPAACESG GTPPSLSHRD ASKELFRVEE DLDQILNLGA EPKPKPQLKP KPPVAAKPVI PRKPAVPPKA GPAAVAVAGQQ KPQEQIQAMD EMDILQYIQD HDTPAQAAPS LF
Source	Mammalian cell
Target Names	HS1BP3
Protein Names	Recommended name: HCLS1-binding protein 3 Alternative name(s): HS1-binding protein 3 Short name= HSP1BP-3
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.