



Recombinant Human Spindlin-4 (SPIN4)

Product Code	CSB-YP704087HU
Abbreviation	SPIN4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q56A73
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSPPTVPPMG VDGVSAYLMK KRHTRKQRR KPTFLTRRNI VGCRIQHGWK EGNEPVEQWK GTVLEQVSVK PTLYIIKYDG KDSVYGLELH RDKRVLALEI LPERVPTPRI DSRLADSLIG KAVEHVFEGE HGTKDEWKGM VLARAPVMDT WFYITYEKDP VLYMYTLLDD YKDGDLRIIP DSNYYFPTAE QEPGEVVDLSL VGKQVEHAKD DGSKRTGIFI HQVVAKPSVY FIKFDDDIHI YVYGLVKTP
Source	Yeast
Target Names	SPIN4
Protein Names	Recommended name: Spindlin-4
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.