



Recombinant Phodopus roborovskii Proto-oncogene c-Fos (FOS)

Product Code	CSB-EP704097PFAG-B
Abbreviation	FOS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q56TN0
Product Type	Recombinant Protein
Immunogen Species	Phodopus roborovskii (Roborovski's desert hamster) (Cricetulus roborovskii)
Purity	≥85% (SDS-PAGE)
Sequence	MMFSGFNTDY EASSSRCSSA SPAGDLSLYY HSPADSFSSM GSPVNAQDFC ADLSVSSANF IPTVTAISTS PDLQWLQPT LVSSVAPSQT RAPHYGVPT PSTGAYSRAG MVKTVSGGRA QSIGRRGKVE QLSPEEEEEKR RIRRERNKMA AAKCRNRRRE LTDTLQAETD QLEDEKSALQ TEIANLLKEK EKLEFILA AH RPACKIPDDL GPFEDMSVAS LDLTGGLPEA ATPESSEAFS LPLLNEPEPK TSLESVKSIS SMELKAEPFD DFLFSASSRP SGSETTARSV PDMDLSGSFY AADWEPLHSS SLGMGPMATE LEPLCTPVVT CTPSCTTYTS SFVFTYPETD SFPSCAAHR KGSSSNPSS DSLSSPTLLA L
Source	E.coli
Target Names	FOS
Protein Names	Recommended name: Proto-oncogene c-Fos Alternative name(s): Cellular oncogene fos
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.