



Recombinant *Corynebacterium ammoniagenes* Riboflavin biosynthesis protein ribF (ribF)

Product Code	CSB-BP704542DWW
Abbreviation	ribF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q59263
Product Type	Recombinant Protein
Immunogen Species	<i>Corynebacterium ammoniagenes</i> (<i>Brevibacterium ammoniagenes</i>)
Purity	>85% (SDS-PAGE)
Sequence	MDIWYGTA AV PKDLNSAVT IGVFDGVHRG HQKLINATVE KAREVGAKAI MVTDFPHVPS VFLPRRAPLG ITTLAERFAL AESFGIDGVL VIDFTRELSG TSPEKYVEFL LEDTLHASHV VVGANFTFGE NAAGTADSLR QICQSRLTVD VIDLLDDEGV RISSTTVREF LSEGDVARAN WALGRHFYVT GPVVRGAGRG GKELGFPTAN QYFHDTVALP ADGVYAGWLT ILPTEAPVSG NMEPEVAYAA AISVGTNPTF GDEQRSVESF VLDRDADLYG HDVKVEFVDH VRAMEKFDSV EQLLEVMKAD VQKTRTLAQ DVQAHKMAPE TYFLQAES
Source	Baculovirus
Target Names	ribF
Protein Names	Recommended name: Riboflavin biosynthesis protein ribF Including the following 2 domains: Riboflavin kinase EC= 2.7.1.26 Alternative name(s): Flavokinase FMN adenylyltransferase EC= 2.7.7.2 Alternative name(s): FAD py
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.