



Recombinant Arabidopsis thaliana B3 domain-containing protein At2g16210 (At2g16210)

Product Code	CSB-MP704727DOA
Abbreviation	At2g16210
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BPT7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MATNAGYLRC KEERKNESFF KVVQSINVSS ENKRALPHDF SRSFTDKELS RKMKIRAQWG NSWEVGIKSN PRFYFMEKSG WEKFVRDNAL GNSELLTFTH KGKMHFTVNI FKLDGKEMMQ PPQSRFFAS SSRIKTEQEE NDIKEEVVVS SNRGQTAAE SKGRKLN LGK RAAKESQSSK RTEKVVRARS DYAGASSSTA AAFTILFKQG YLVFLRIPNS VSKDQVPDEK TVFKIHHPNG KKSWNVVYLE RFGAFSGGWR RVVKEYPLAV GDTCKFTFIK PKELLLVSK P
Source	Mammalian cell
Target Names	At2g16210
Protein Names	Recommended name: B3 domain-containing protein At2g16210
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.