



Recombinant Mouse Endoribonuclease ZC3H12A (Zc3h12a)

Product Code	CSB-MP704737MO
Abbreviation	Zc3h12a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5D1E7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSDPCGTKPVQESNPTMSLWSLEDRHSSQGRPQPDQDPVAKEAPTSELQMK VDFFRKLG SSSEIHSV LQKLG VQADTNTVLGELVKHGSATERECQALTAPSPQPPLVPRGG STPKPST LEPSLPEEDREGSDLRPVVIDGNSVAMSHGNKEVFSCRGILLAVNWFLERGH TITVFVP SWRKEQPRPDVPITDQHILRELEKKKILVFTPSRRVGGKRVVCYDDRFIVKLAF ESDGVV VSNDTYRDLQGERQEWKRFIEERLLMYSFVNDKFMPPDDPLGRHGPSLDNFL RKKPLPSE HRKQPCPYGKKCTYGIKCRFFHPERPSRQRSVADEL RANALLSPRTPVKD KSSQRPS ASQSSSVSLEAEPGSLDGKKGARSSPGPHREGSPQTCAPAGRSLPVSGGSF GPTWLAH TQDSLPTYTSQECLDSGIGSLESQMSELWGVRGGSPGESGPTRGPYAGYHSY GSKVPAAPS FSPFRPAMGAGHFSVPTDYVPPPPTYP SREYWSEPYLPPPTPVLQEPQRPS PGAGGGPW GRVGD LAKERAGVYTKLCGVFPPHLVEAVMRRFPQLLDPQQLAAEILSYKSQH LSE
Source	Mammalian cell
Target Names	Zc3h12a
Protein Names	Recommended name: Ribonuclease ZC3H12A EC= 3.1.-.-Alternative name(s): MCP-induced protein 1 Zinc finger CCCH domain-containing protein 12A
Expression Region	1-596
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.



Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.