



Recombinant Bovine Complement C1q subcomponent subunit A (C1QA)

Product Code	CSB-EP704831BO
Abbreviation	C1QA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5E9E3
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	EDVCRAPD GTHGSAGIPG RPGRPGLKGE RGEPEGAPAIQ TGIRGLKGDQ GDPGPPGNPG RMGYPGPSGP MGPAGLPGLK GTKGSPGNIK DQPRPAFSAV GPNSVSRDNV VVFGKVITNQ ENVYQNNTGR FRCSVPGYYY FTFQVVSNWD ICL SIRSSRR DQIQPLGFCD FNSKGGFFQVV SGGTVLHLQQ GDQVWIEKDP SKGRIYHGSE ADSIFSGFLI FPSPA
Source	E.coli
Target Names	C1QA
Protein Names	Recommended name: Complement C1q subcomponent subunit A
Expression Region	23-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a major constituent of the human complement subcomponent C1q. C1q associates with C1r and C1s in order to yield the first component of the serum complement system. Deficiency of C1q has been associated with lupus erythematosus and glomerulonephritis. C1q is composed of 18 polypeptide chains: six A-chains, six B-chains, and six C-chains. Each chain contains a collagen-like region located near the N terminus and a C-terminal globular region. The A-, B-, and C-chains are arranged in the order A-C-B on chromosome 1. This gene encodes the A-chain polypeptide of human complement subcomponent C1q.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.