



Recombinant Bovine Nitrogen permease regulator 2-like protein (NPRL2)

Product Code	CSB-EP704841BO-B
Abbreviation	NPRL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5E9U9
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MGSSCRIECI FFSEFHPTLG PKITYQVPED FISRELFDTV QVYIITKPEL QNKLITVTAM EKKLIGCPVC IEHKKYSRNA LLFNLGFVCD AQAKTCALEP IVKLAGYLT TLELESSFVS TEESKQQLVP IMTILLEELN ASGRCTLPID ESNTIHLKVI EQRPDPPVAQ EYDVPVFTKD KEDFFNSQWD LTTQQILPYI DGFRHVQKIS AEADVLENLV RIAIQNLLYY GVVTLVSILQ YSNVYCPTPK VQDLVDDKSL QEACLSYVTK EGHKRASLRD VFQLYCSLSP GTTVRDLIGR HPQQQLQHVDE RKLIQFGLMK NLIRRLQKYP VRVSRDERSH PARLYTGCHS YDEICCKTGM SYQELDERLE NDPNIIICWK
Source	E.coli
Target Names	NPRL2
Protein Names	Recommended name: Nitrogen permease regulator 2-like protein Short name= NPR2-like protein Alternative name(s): Tumor suppressor candidate 4
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.