




Recombinant Mouse Interferon regulatory factor 3 (Irf3)

Product Code	CSB-YP011818MO
Relevance	Key transcriptional regulator of type I interferon (IFN)-dependent immune responses which plays a critical role in the innate immune response against DNA and RNA viruses. Regulates the transcription of type I IFN genes (IFN-alpha and IFN-beta) and IFN-stimulated genes (ISG) by binding to an interferon-stimulated response element (ISRE) in their promoters. Acts as a more potent activator of the IFN-beta (IFNB) gene than the IFN-alpha (IFNA) gene and plays a critical role in both the early and late phases of the IFNA/B gene induction. Found in an inactive form in the cytoplasm of uninfected cells and following viral infection, double-stranded RNA (dsRNA), or toll-like receptor (TLR) signaling, is phosphorylated by IKBKE and TBK1 kinases. This induces a conformational change, leading to its dimerization and nuclear localization and association with CREB binding protein (CREBBP) to form dsRNA-activated factor 1 (DRAF1), a complex which activates the transcription of the type I IFN and ISG genes. Can activate distinct gene expression programs in macrophages and can induce significant apoptosis in primary macrophages.
Abbreviation	Recombinant Mouse Irf3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P70671
Alias	Short name:IRF-3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	METPKPRILPWLVSQDLGQLEGVAWLDESRTFRIPWKHGLRQDAQMADFG IFQAWAEASGAYTPGKDKPDVSTWKRNFERSALNRKEVLRRLAADNSKDPYDPH KVYEFVTPGARDFVHLGASPDTNGKSSLPHSQENLPKLFDFGLILGPLKDEGSS DLAIVSDPSQQLPSPNVNNFLNPAPQENPLKQLLAAEEQWFEVTAFYRGRQV QQTLFCPGGLRLVGSTADMTLPWQPVTLPDPEGFLTDKLVKEYVGGVVKGLG NGLALWQAGQCLWAQRLGHSHAFWALGELLPSGRGPDGEVHKDKDGAV FDLRPFVADLIAFMEGSGHSPRYTLWFCMGEMWPQDQPWVKRLVMVKVPT CLKELLEMAREGGASSLKTVDLHISNSQPISLTSQYKAYLQDLVEDMDFQAT GNI
Research Area	Immunology
Source	Yeast
Target Names	Irf3



Protein Names	Recommended name: Interferon regulatory factor 3 Short name= IRF-3
Expression Region	1-419aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	50.9kDa
Protein Length	Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.