



Recombinant *Oryza sativa* subsp. *japonica* Protein kinase and PP2C-like domain-containing protein (Os01g0541900, LOC_Os01g36080), partial

Product Code	CSB-BP705502OFG
Abbreviation	Os01g0541900, LOC_Os01g36080
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JJY4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Source	Baculovirus
Target Names	Os01g0541900
Protein Names	Recommended name: Protein kinase and PP2C-like domain-containing protein Including the following 2 domains: Probable serine/threonine-protein kinase Os01g0541900 EC= 2.7.11.1 Probable protein phosphatase 2C 4 Short name= OsPP
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.