



# Recombinant Human PCI domain-containing protein 2 (PCID2)

<b>Product Code</b>	CSB-EP705514HU
<b>Abbreviation</b>	PCID2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5JVF3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AHITINQYL QQVYEAIISR DGASCAELVS FKHPHVNPR LQMASPEEKC QQVLEPPYDE MFAAHLRCTY AVGNHDFIEA YKCQTVIVQS FLRAFQAHKE ENWALPVMYA VALDLRVFAN NADQQLVKKG KSKVGDMLEK AAELLMSCFR VCASDTRAGI EDSKKWGMLF LVNQLFKIYF KINKLHLCKP LIRAIDSSNL KDDYSTAQRV TYKYYVGRKA MFDSDFKQAE EYLSFAFEHC HRSSQKNKRM ILIYLLPVKM LLGHMPTVEL LKKYHLMQFA EVTRAVSEGN LLLLHEALAK HEAFFIRCGI FLILEKLIKII TYRNLFKKVY LLLKTHQLSL DAFLVALKFM QVEDVDIDEV QCILANLIYM GHVKGYISHQ HQKLVVSKQN PFPPLSTVC
<b>Source</b>	E.coli
<b>Target Names</b>	PCID2
<b>Protein Names</b>	Recommended name: PCI domain-containing protein 2 Alternative name(s): CSN12-like protein
<b>Expression Region</b>	2-399
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.