



Recombinant *Oryza sativa* subsp. *japonica* Expansin-A17 (EXPA17)

Product Code	CSB-BP706017OFG
Abbreviation	EXPA17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4PR49
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	AQVVTAGFT TDLYWQQQPA PGAVTPYKTS DWHDGSATFY GDPSGMGDDF GGACGYVSND IVSLYSTKTA ALSTPLFADG NGCGQCYELR CVKSPWCNPG SPSVVITGTN LCPPNWYLPN DDGGWCNPPR HHFDMAPPSF LKLAQRVAGI VPVQYRRVPC QRTGGVRFCL QGNHYWLLLY VMNVGGAGDV SSLSVKTSGG GGAWIQAAHN WGITYQVFAA LDNSDGLTVK LTTYSTPQQT IIVSDAISPW WITGLCYQGS NNFY
Source	Baculovirus
Target Names	EXPA17
Protein Names	Recommended name: Expansin-A17 Alternative name(s): Alpha-expansin-17 OsEXP17 OsEXPA17 OsaEXPa1.31
Expression Region	22-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.