



# Recombinant *Macaca fascicularis* Glutamine synthetase (GLUL)

<b>Product Code</b>	CSB-EP706159MOV-B
<b>Abbreviation</b>	GLUL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q4R7U3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus</i> monkey)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TTSASSHLN KGIKQVYMSL PQGEKVQAMY IWIDGTGEGL RCKTRTL DSE PKC VEELPEW NFDGSSTLQS EGSNSDMYLV PAAMFRDPFR KDPNKLVLCE SFQVQFEGPA ETNLRHTCKR IMDMVSNQHP WFGMEQEYTL MGTDGHPFGW PSNGFPGPQG PYYCGVGADR AYGRDIVEAH YRACLYAGVK IAGTNAEVMP AQWEFQIGPC EGISMGDHLW VARFILHRVC EDFGVIATFD PKPIPGNWNG AGCHTNFSTK AMREENGLKY IEEAIEKLSK RHQYHIRAYD PKGGLDNARR LTGFHETSNI NDFSAGVANR SASIRIPRTV GQEKKGYFED RRPSANCDPF SVTEALIRTC LLNETGDEPF QYKN
<b>Source</b>	<i>E. coli</i>
<b>Target Names</b>	GLUL
<b>Protein Names</b>	Recommended name: Glutamine synthetase Short name= GS EC= 6.3.1.2 Alternative name(s): Glutamate decarboxylase EC= 4.1.1.15 Glutamate-- ammonia ligase
<b>Expression Region</b>	2-374
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.