



Recombinant Danio rerio Inositol oxygenase (miox)

Product Code	CSB-BP706312DIL
Abbreviation	miox
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4V8T0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MGPDP SLAYR PECHEKDKTE FRNFENGDLF DRVFN TYKLM HTHQTLD FVK QKHQVWSNCS HFSLSMMDSI DSLDELVD ESDPDVDFPNSF HAFQTAEGIR REHPDKDWFQ LVGLIHDVGK VMALYSEPQW AVVGDTYPVG CKFQNSIVFR NSTFEGNPDG KNPAPNTEFG IYEPQCGLDK VLMSWGHDEY LYRVMKFENKC TIPEEGLYMI RFHSFY PWHS NGDYMHL CNE KDQQMLPWWK EFNKFDLYTK STELPDVERL KPYYQSLIDK YCPGVLQW
Source	Baculovirus
Target Names	miox
Protein Names	Recommended name: Inositol oxygenase EC= 1.13.99.1 Alternative name(s): Myo-inositol oxygenase Short name= MI oxygenase
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.