



Recombinant Guanine nucleotide-binding protein alpha-14 subunit (gpa-14)

Product Code	CSB-MP706326CXX
Abbreviation	gpa-14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4VT38
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	>85% (SDS-PAGE)
Sequence	MFSCFNNLGL DYCQCMHGP EGCMVPSRQG DGGELYAHSE ELEAKLRGLA KKELEIEKS LENDKKTYS HIKILILGGP SSGKSTIFKQ MQIIHSNGFK TEQELIQYRG LIDTNRQTY RQLVSGARVV GISLESLESL VHDINKVYAP MAADEFSIRT IPDVVEPLTA FWNSREIQEV YKRRYEFELL DSTKYYLENL NRISKSDYLP NEEDIVHSRK ATVSINSIVF QYTGVSLLMV DVGGQRSERK KWLHLFDDAK VVIFVIDLTG YAKKSEESRT ELSRFPNFFN EIGNDAFDMK VALKIFNDVA GSHALANAVF LFFNKVDLF KELLPQVSLQ PCFSKFAEEN SYDNTSKFIC DKFIRAAKPK KSVFPHFTTA TNTENIKMVF RACMESVFKA NSKATGLS
Source	Mammalian cell
Target Names	gpa-14
Protein Names	Recommended name: Guanine nucleotide-binding protein alpha-14 subunit
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.