



Recombinant Danio rerio Cyclin-related protein FAM58A (fam58a)

Product Code	CSB-BP706502DIL
Abbreviation	fam58a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q503D6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MSAHASTSAA ANEARGRRSA EDVEYSKTHF RVCRFITETG VKLGMRSVPM ATACVLYHRF FQSASLQIYE PYLVAMSAIH LAGKVEEQHL RTRDIINVCH RYFHPDSEPL ELNGKFWELR DSIVQCELLI LRQLNFQVTF EHPHKYLLHY LLSVRLLNR HAWSRTPIAE TALAVLKDSY HGSVCVRHRP QHLALTALYL ALQTYGVQLP RGELEWWQVV CADITKAQIE TIMSELLQLY DMEAKCT
Source	Baculovirus
Target Names	ccnq
Protein Names	Recommended name: Cyclin-related protein FAM58A
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.