



Recombinant Human Pleckstrin homology domain-containing family O member 1 (PLEKHO1)

Product Code	CSB-EP706633HU-B
Abbreviation	PLEKHO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53GL0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMKKNNSAKR GPQDGNQQPA PPEKVGWVRK FCGKGIFREI WKNRYVVLKG DQLYISEKEV KDEKNIQEVF DLSDYЕКCEE LRKSKRSRKK NHSKFTLAHS KQPGNTAPNL IFLAVSPEEK ESWINALNSA ITRAKNRILD EVTVEEDSYL AHPTRDRAKI QHSRRPPTRG HLMVAVASTST SDGMLTLDLI QEEDPSPEEP TSCAESFRVD LDKSVAQLAG SRRRADSDRI QPSADRASSL SRPWEKTDKG ATYTPQAPKK LTPTEKGRCA SLEEILSQRD AASARTLQLR AEEPPTPALP NPGQLSRIQD LVARKLEETQ ELLAEVQGLG DGKRKAKDPP RSPDSESEQ LLETERLLG EASSNWSQAK RVLQEVREL R DLYRQMDLQT PDSHLRQTPP HSQYRKSLM
Source	E.coli
Target Names	PLEKHO1
Protein Names	Recommended name: Pleckstrin homology domain-containing family O member 1 Short name= PH domain-containing family O member 1 Alternative name(s): C-Jun-binding protein Short name= JBP Casein kinase 2-interacting protein 1 Sho
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.