



Recombinant *Synechococcus elongatus* Amidophosphoribosyltransferase (purF)

Product Code	CSB-YP706926FPY
Abbreviation	purF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55038
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	CGVF ALYAPGEEVA RMAYFGLYAL QHRGQESAGI AVFEGDRVML HKDMGLVSQV FDPEILQQLQ GSLAVGHTRY STTGSSRIAN AQPALLETRL GPVALAHNGN LVNTVELRQE LLAKNHELTT TTDSELIAFA IMEAVAEGQD WRGAIESACR RSQGAFSLTI GTPEALYGTR DPNGIRPLVL GTLESNGQTR YVLSSETCGL DIIGADYVRD IAPGEMVRIT DAGLESWTWA EAPQPKLCVF EMIIYFARPDS LFHGESLYSY RRRIGQRLAK EAPADVDLVL GVPDSGIPAA IGFSEALGIP YAEGLIKRY VGRTFIQPTQ SMRETGIRMK LNPLKDVLAG KRIAIIDDSI VRGTTSRKLV KALRDAGATE VHMRISSPPV THPCFYGIDT DTQDQLIAAT RSVSEITEQI GVDSLAYLTE QGMLEATRES IGNFCTACFN GRYPPIAIPPEE IKRSKLMLET VTA
Source	Yeast
Target Names	purF
Protein Names	Recommended name: Amidophosphoribosyltransferase Short name= ATase EC= 2.4.2.14 Alternative name(s): Glutamine phosphoribosylpyrophosphate amidotransferase Short name= GPATase
Expression Region	27-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.