



Recombinant Dictyostelium discoideum GTP-binding protein Rheb homolog (rheb)

Product Code	CSB-MP706931DKK
Abbreviation	rheb
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q550Q4
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MAPQKHRKIC VMGSRAVGKS TITMQFVESH CPDSYHPTIE NTYQKIIRHQ GQEYSVEIID TAGQDEYSIL QKQYSIGIHG YILVYSVTSV SSLEVIKVLN DKILSSLGAE KIPRVLVGNK SLDNERNIS RDQGIQLANE WECAFVECSG KNNENVEEVF KQILNEVNKG STGPEPPQKE GC
Source	Mammalian cell
Target Names	rheb
Protein Names	Recommended name: GTP-binding protein Rheb homolog
Expression Region	1-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.