



Recombinant *Synechococcus elongatus* Lycopene beta cyclase (crtL)

Product Code	CSB-BP706938FPY
Abbreviation	crtL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55276
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	MFDALVIGSG PAGLAIAAEL AQRGLKVQGL SPVDPFHPWE NTYGIWGPPEL DSLGLEHLFG HRWSNCVSYF GEAPVQHQYN YGLFDRAQLQ QHWLRQCEQG GLQWQLGKAA AIAHDSHHSC VTTAAGQELQ ARLVVDTTGH QAAFIQRPHS DAIAYQAAYG IIGQFSQPPI EPHQFVLMDY RSDHLSPEER QLPPTFLYAM DLGNDVYFVE ETSLAACPAI PYDRLKQRLY QRLATRGVTV QVIQHEEYCL FPMNLPLPDL TQSVVGFSGA ASMVHPASGY MVGALLRRAP DLANAIAAGL NASSSLTAE LATQAWRWGLW PTEKIRKHYI YQFGLEKLMR FSEAQLNHHF QTFGLPKEQ WYGFLTNTLS LPELIQAMLR LFAQAPNDVR WGLMEQQGRE LQLFWQAIAA R
Source	Baculovirus
Target Names	crtL
Protein Names	Recommended name: Lycopene beta cyclase EC= 5.5.1.19
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.