



Recombinant Dictyostelium discoideum WASH complex subunit CCDC53 homolog (DDB_G0267948)

Product Code	CSB-MP707041DKK
Abbreviation	DDB_G0267948
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55FU3
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	≥85% (SDS-PAGE)
Sequence	MQVYVDYRDV NSISNKKTLA YVNNFVINTT QFLNRFNYLC EQKLCDISNH LQKLEITMSL LEIKLGSIDG LESLPSSANV PVHNTQPTTQ QNTPPPPTAA ATTSSNIPPP PPPPPPPMTG VPPPPPPPPP PPIKSNIPP PPATQQNTQE SDDDDDEDNNN NDKNNNNRDD EDDDDDETSSQ APPPPPMPTL PSSQSLMDNL AAVLAMKKQQ RNPQPPQPQP QSPSPQPPPP PTTTSSIPVP PPPFSGANNN SDDDDDDDEG PLSLIPAPPS MMVMSGNNNS NSNSYSNNNN NKDDDDDDDD DDDDDDNSIS IPKPPQMDDI MIPAPPTDLF IPAPPSIRNY DDADDDDDDD DESW
Source	Mammalian cell
Target Names	washc3
Protein Names	Recommended name: WASH complex subunit CCDC53 homolog
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.