



Recombinant Ureaplasma parvum serovar 3 Urease accessory protein UreD (ureD)

Product Code	CSB-MP707074UAO
Abbreviation	ureD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q56562
Product Type	Recombinant Protein
Immunogen Species	Ureaplasma parvum serovar 3 (strain ATCC 700970)
Purity	>85% (SDS-PAGE)
Sequence	MILNKEKIKN YAAYLYIKVA YDQAHSKMAH TVYFTNFYRS SKPLFLDEED PINPCFQTIS MGGGYVSGEI YRSDFEINDD ARCIITTQSS AKAYKTVDGK TSEQHTNITL GKNSILEYIS DNVIVYEDGK FAQFNNFKMD SSATLIYTEC FGPGWSPHGS AYQYEKMYLN TKIYYDDKLV LFDNLKFQPR KNDESAFGIM DGYHYCGTMI VINQQVIEDD VIKIRDLVKE KYPDMDMIFG VSRMDIPGLG LRVLANTYYH VEKINAVAHD YFRRKLFNKK PLILRKP
Source	Mammalian cell
Target Names	ureD
Protein Names	Recommended name: Urease accessory protein UreD
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.