



Recombinant Danio rerio Protein odd-skipped-related 2 (osr2)

Product Code	CSB-BP707085DIL
Abbreviation	osr2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q567J8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MGSKTLPAPV PLHPAPQLNY SLLRTLNAFP VAVDQLYGLS ALHTVQMNRW TVGFPQLQGL ADPRFPGALP FPAAAAHLLP HKHPVHRKDR PRFDFANLAV AATQEDPPVT GQSRLSPERR PARGRLPAKS KKEFICRFCG RHFTKSYNLL IHERHTTDER PYTCDICHKA FRRQDHLRDH RYIHSKEKPF KCQECGKGFC QSRTLAVHKT LHLQTSSLISA AAAAGRCFGA TATCGATV
Source	Baculovirus
Target Names	osr2
Protein Names	Recommended name: Protein odd-skipped-related 2 Short name= zOsr2
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.