



Recombinant Danio rerio N-alpha-acetyltransferase 40 (naa40)

Product Code	CSB-EP707090DIL-B
Abbreviation	naa40
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q568K5
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	GRKSNRAKE KKQRRLEERA AMDAVCAKVD AANKLEDPLS AMPVFKKYDR NGLNLQIECK RVTALSPDTV EWAYELTRAN MQTLYEQSEW GWKEREKREE MKDERAWYLL ARDADSTPLA FSHFRFDVEC GDEVLYCYEV QLESKVRKRG LGKFLIQILQ LIANSTQMKK VMLTVFKHNNH GAYQFFREAL QFEIDETSPS MSGCCGEDCS YEILSRRTKY GEVSGHAHGG HCGGCCH
Source	E.coli
Target Names	naa40
Protein Names	Recommended name: N-alpha-acetyltransferase 40 EC= 2.3.1.- Alternative name(s): N-acetyltransferase 11 NatD catalytic subunit
Expression Region	2-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.