



Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At5g35995 (At5g35995)

Product Code	CSB-BP707115DOA
Abbreviation	At5g35995
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q56W59
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MGSRDFISSL PDEVLGKKIL SLLPTKLVVS TSVLSKRWRN LFHFVDSFDL EDSTPIRNAD SFSDFMKDTV ALLSNCPIKR LTLNSHYEKS SVNRWIRSAL KRGCLELNLQ SQYDHYLDIG IFFRNNTLVK LTLSSYRTFL RGNVPPEGRV FFPALKTSL GAVVAKPALY NWLISGCPVL EELFILDVGS GDDQPTWTRS VVSASIKRLT IVFHFPHTY PYEDDVEIKT PNLEFLDYSA LRSDGSDVDY LDSLAEARLD LRLWELTNTF DFGEVTNLVS AIRNVKTLHL SSSSLEAFYY RCYTMPVFDK LIHLSIESDK ENGWQALPRL LLKSPNLQTL AIKGLLHKVT YRCGNACACT SKPKKRYLYK RYEDDSPTSE VEGRCCLSTC RVKVLISGY GGSLREVKQM EHFLRKLKYL ETVKIGVEED NDSEHLRTEL MTLDIASSKC NIQFI
Source	Baculovirus
Target Names	At5g35995
Protein Names	Recommended name: F-box/LRR-repeat protein At5g35995
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.