



# Recombinant Arabidopsis thaliana Patellin-3 (PATL3)

<b>Product Code</b>	CSB-EP707121DOA-B
<b>Abbreviation</b>	PATL3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q56Z59
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAEPTTTTTL VTPEKLPSPS LTPSEVSEST QDALPTETET LEKVTETNPP ETADTTTKPE EETA AEHPP TVTETETAST EKQEVKDEAS QKEVAEEKKS MIPQNLGSFK EESSKLS DLS NSEKKS LDEL KHLVREALDN HQFTNTPEEV KIWGIPLLED DRSDVLLKF LRAREFKVKD SFAMLKNTIK WRKEFKIDEL VEEDLVDDLD KVVFMHGHDR EGHPVCYNVY GEFQNKELYN KTFSDDEEKRK HFLRTRIQFL ERSIRKLDFS SGGVSTIFQV NDMKNSPGLG KKELRSATKQ AVELLQDNYP EFVFKQAFIN VPWWYLVFYT VIGPFMTPRS KSKLVFAGPS RSAETLFKYI SPEQVPVQYG GLSVDPCDCN PDFSLED SAS EITVKPGTKQ TVEIIYEKC ELVWEIRVTG WEVSYKAEFV PEEKDAYTVV IQKPRKMRPS DEPVLTHSFK VNELGKVLLT VDNPTSKKKK LVYRFNVKPL
<b>Source</b>	E.coli
<b>Target Names</b>	PATL3
<b>Protein Names</b>	Recommended name: Patellin-3
<b>Expression Region</b>	1-490
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.