



Recombinant Arabidopsis thaliana Patellin-3 (PATL3)

Product Code	CSB-MP707121DOA
Abbreviation	PATL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q56Z59
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAEPTTTTTL VTPEKLPSPS LTPSEVSEST QDALPTETET LEKVTETNPP ETADTTTKPE EETAAEHHPP TVTETETAST EKQEVKDEAS QKEVAEEKKS MIPQNLGSFK EESSKLSDLN NSEKKSLEDEL KHLVREALDN HQFTNTPEEV KIWGIPLLED DRSDVLLKF LRAREFKVKD SFAMLKNTIK WRKEFKIDEL VEEDLVDDLD KVVFMHGHDR EGHPVCYNVY GEFQNKELYN KTFSDDEEKRK HFLRTRIQFL ERSIRKLDLFS SGGVSTIFQV NDMKNSPGLG KKELRSATKQ AVELLQDNYP EFVFKQAFIN VPWWYLVFYT VIGPFMTPRS KSKLVFAGPS RSAETLFKYI SPEQVPVQYG GLSVDPCDCN PDFSLEDSAS EITVKPGTKQ TVEIIYEKC ELVWEIRVTG WEVSYKAEFV PEEKDAYTVV IQKPRKMRPS DEPVLTHSFK VNELGKVLTL VDNPTS KKKK LVYRFNVKPL
Source	Mammalian cell
Target Names	PATL3
Protein Names	Recommended name: Patellin-3
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.