



Recombinant Candida albicans RuvB-like helicase 2 (RVB2)

Product Code	CSB-MP707592CZD
Abbreviation	RVB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5AGZ9
Product Type	Recombinant Protein
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MATNTITTTK VTTTTTDVTG LSLIAAHSHI SGLGLDDNLQ PKENAQGMVG QLSARKAAGV ILKMVEAGKI AGRAVLIAGP PSTGKTAIAM GLSQSLGNQV PFTALAASEV FSLELSKTEA LTQAFRKSIG IKIKEETEI EGEVVEIQID RTITGGHKQG KLTIKTTDME TIYELGNKMI EGLTKEKVLA GDVISIDKAS GKITKLGRSF TRARDYDAMG PETKFVQCPE GELQKRKEVV HTISLHEIDV INSRQQGFLA LFSGDTGEIR PEVRDQINTK VAEWKKEGKA EIVPGVLFID EVHMLDIECF SFINRALEDE FSPIVIMATN RGVSRIKGTD YKSPHGMPPMD LLDRSITIHT TSYTADEIRT ILSIRATEEE VELSGDALAL LTKIGQETSL RYAANLISVS QQIALKKKNN TVDLQDIKRA YMLFLDSDRS VQYLEENADQ YIDDYGRVTI GQUESTDGSTQ PQAKQQEVAQ PEATQPQSQP EDDKMETD
Source	Mammalian cell
Target Names	RVB2
Protein Names	Recommended name: RuvB-like helicase 2 EC= 3.6.4.12
Expression Region	1-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.