



Recombinant Danio rerio Uncharacterized protein C6orf106 homolog (zgc:101577)

Product Code	CSB-EP707686DIL
Abbreviation	zgc:101577
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BL31
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEGMDIDLQ ELMQKFSCMG TTDKDLVSE FQRLGFLN PAGCAFFLDM TNWNLQAAIG AYYDFESPNI NTPSMSFVED VTIGEGESVP PDTPFTKTWR IQNTGTESWP PGVCLKYVGG DQFGHVMVM VRSLDPQEIS DVSVMRSPA VPGMYQGQWR MCTATGLFYG DVIWVLSVE EGLLGVTQQ LSSFKTEFNT QPHRSLEGDY NPFASPQKNK QDTNEDHLKD PGGPWEAPLD SIQQDQNGLN HSSVNITPNG LQNNLSVVTY SQGINGPFPF GQS
Source	E.coli
Target Names	zgc:101577
Protein Names	Recommended name: Uncharacterized protein C6orf106 homolog
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.